

Date d'édition : 25.05.2026

Ref : 662088

Lot de produits chimiques colorants



Liste de livraison

- 100 ml
Methylene blue after• LÖFFLER, for staining of fungi, bacteria and nuclei
 - 100 ml
Carbolic gentian violet solution for GRAM staining
 - 100 ml
Iodine-potassium iodide solution after LUGOL
 - 100 ml
Azure-eosine-methylene-blue solution after GIEMSA for staining of blood and protozoa•••
 - 25 g
Eosine, yellowish for unspecific staining•••••
 - 1 l
2-propanol, for preparation of alcohol dilution series
 - 250 ml
Xylene
 - 50 ml
Embedding medium
- Import texte : janvier 2015

Catégories / Arborescence

Sciences > Biologie > Produits > Microscopy > Accessories for microscopy